

How Effective is your Temperature Monitoring System?

Need an accurate and reliable temperature recording system notifying you when your drug or sample storage temperature is out of range?

Introducing the new AS300SA GSM Alarm Monitoring system:



- Monitoring temperature ranges between -200°C to +100°C
- Simple Installation
- Setup or change settings from a mobile phone in any country
- 4000 reading memory
- Dial/ Message up to 5 telephone numbers (mobiles/ landlines) when temperature is out of range and alarm occurs.
- Default user can dial into system for readings. Download details include: alarm status, signal strength, setup values, and readings
- Battery Backup
- Calibration Function
- Test alarm Function
- Make an informed decision on what action to take based on instant readings, 24 hours a day 7 days a week.

- 24 Hour Support
- Calibration and Repair Service
- End of Life Recycling
- New Technology
- Hire and Purchase options
- Global Procurement
- Lowest Total Cost Solution for Your Clinical Study

Your Global Partner of choice for Clinical Study Equipment



tel. +44 (0) 1695 721980
fax. +44 (0) 1695 721990

sales@mesm.co.uk www.mesm.co.uk